




21 CFR Part 11 Compliance Readiness

User administration
(Local oder active directory authentication)




Comprehensive activity logging (by object, by user)

The screenshot shows an "Audit trail of Cannabis 3D" window with a table of activity logs. The table has columns for Date, Person, Location, and Change.

Date	Person	Location	Change
08.06.2006 14:46:00	ks	lab007	Imported from 'C:\Documents and Settings\ks.LABCOGNITION\My Documents\My panorama
08.06.2006 14:46:24	ks	lab007	Data name changed from 'Testmessung; Path 1 (0,00mm)' to 'Cannabis 3D'
26.02.2008 12:14:15	tester	QA	Derivative: Order='1', Algorithm='Savitzky-Golay', Smoothing Window='Nine Points'
26.02.2008 12:14:30	tester	QA	Normalization: Maximum='1', Minimum='0', Normalize Minimum='Yes'

Electronic signature & audit trail for data objects

Data and functional access permission schemes



Release and unrelease objects in keeping with organizational data control strategies